

شركة المجهر المتكامل التجارية

AL MEJHAR AL MOTAKAMEL COMMERCIAL CO.

Benchtop Centrifuge
1790xg Centrifugal Force
20mlx6 Rotor Capacity



شركة المجهر المتكامل التجارية

AL MEJHAR AL MOTAKAMEL COMMERCIAL CO.

2 Knobs Operation

Control time and speed

Non slip rubber pads for stable operation



 Centrifuge



800



800D



800-1



800 digital



Product Usage

The frame of this instrument is made of metal .Its model is beautiful, and it has the advantages of small volume, low weight,large capacity,low noise,high efficiency and so on.

It can be used in hospitals and biochemical labs for qualitative analysis of serum, urea and plasma.

Features

- Separating the liquid from the solid particles or the liquid mixture of different density;
- Equipped with speed and timer control;
- Widely used in the laboratory, medical, chemistry and so on.

Model	800	800D/800-1	800digital	80-1	80-2
Voltage	AC 220V/110V 50Hz/60Hz	AC 220V/110V 50Hz/60Hz	AC 220V/110V 50Hz/60Hz	AC 220V/110V 50Hz/60Hz	AC 220V/110V 50Hz/60Hz
Max speed	4000rpm	4000rpm	4000rpm	4000rpm	4000rpm
Speed indication	Stepless speed adjustable	Stepless speed adjustable	Stepless speed adjustable, Digital display	Stepless speed adjustable	Stepless speed adjustable
Capacity	20mlx6 10mlx6 (angle rotor)	20mlx6 10mlx6 (angle rotor)	20mlx6 10mlx6 (angle rotor)	20mlx6 10mlx8 (angle rotor)	20mlx12 10mlx8 (angle rotor)
Max rcf	1790xg	1790xg	1790xg	1790xg	2325xg
Timer	----	0-60min	0-60min	0-120min	0-120min
Dimension	Φ 250x220mm	Φ 250x220mm	Φ 250x220mm	Φ 245x250x220mm	Φ 280x305x270mm
Package size	30x30x27cm	30x30x27cm	30x30x27cm	32x30x28cm	37x34x33cm
Weight	4kg	4kg	4kg	6kg	8.5kg



Low speed Laboratory centrifuge

800 electric centrifuge is widely used in hospitals, chemical and Biological Chemistry Laboratory for serum, plasma, immune qualitative analysis experiment. Has beautiful shape, large capacity, small volume, complete function. Its performance is stable, the speed is adjustable and can automatically adjust the balance, low temperature, high efficiency and wide applicability and advantages.



Type No.	800	800D/800-1	800-3 digital display
Speed range	0~4000 (rpm)		
Capacity (ml)	20 (ml) ×6		
Time Range (min)	/	0~120 (min)	
Maximum centrifugal force (xg)	1790 (xg)		
Voltage	220V 50HZ		110V 60HZ
Package size	31*31*28CM		
Material	Metal		
Weight	4KG		

شركة المجهر المتكامل التجارية

AL MEJHAR AL MOTAKAMEL COMMERCIAL CO.



0-4000RPM Speed

Used for separating suspensions and liquids