



# AL MEJHAR AL MOTAKAMEL CO.

## 80-2B China Medical Equipment Laboratory Plasma Analysis Centrifuge



### Technical specification:

product name:	low speed centrifuge in medical
model:	80-2B
speed range:	start up-4000rpm, non-step adjusted speed
capacity:	20ml x12
Max RCF:	2325 x g
time range:	0-30min

### GENERAL

The frame of this instrument is made of high strength plastic. Its model is beautiful, and it has the advantages of small cubage, low weight, large capacity, low noise, high efficiency and so on. It can be used in hospitals and biochemical labs for qualitative analysis of serum, urea and plasma.

### PECULIARITY

- Rotor imbalance protection device
- Speed choice in advance & acceleration accurate gradually
- Door switch

